

**AMENDMENT NO. 1 MAY 2002  
TO  
IS 4028 : 1992 BEESWAX, BLEACHED FOR COSMETIC  
INDUSTRY — SPECIFICATION  
( Third Revision )**

[ Page 1, Table 1 Sl No. (ii) ] — Substitute the following for the existing:

| (1) | (2)        | (3)  | (4) |
|-----|------------|------|-----|
| (i) | Acid value | 5-15 | A-3 |

[ Page 1, Table 1, Sl No. (iii), col 3 ] Substitute '87-104' for '85-95'.

[ Page 1, Table 1, Sl No. (viii)] — Insert the following after Sl No. (viii):

| (1) | (2)                                       | (3)                   | (4)  |
|-----|-------------------------------------------|-----------------------|------|
| ix) | Glycerin and other<br>polyhydric alcohols | To pass the test<br>✓ | A-10 |

( Page 3, clause A-4.2.4 ) — Substitute 'Indicator' for 'Indication'.

( Page 5, clause A-8.2, line 3 ) — Substitute '5.0 g' for '50 g'.

( Page 6, clause A-9.2 ) — Insert the following new clause after A-9.2:

**A-10 TEST FOR GLYCERIN AND OTHER POLYHYDRIC ALCOHOLS**

**A-10.1 Reagents**

**A-10.1.1 Alcoholic Potassium Hydroxide Solution** — approximately 0.5 N, prepared by dissolving potassium hydroxide in rectified spirit.

**A-10.1.2 Sulphuric Acid** — 1 M.

**A-10.1.3 Sodium Periodate** — 1.07 percent m/v.

**A-10.1.4 Glycerin** — 0.001 percent m/v.

**A-10.1.5 Fuchsin Solution**

**Amend No. 1 to IS 4028 : 1992**

**A-10.2 Procedure**

Take 0.2 g of the sample and add 10 ml of alcoholic potassium hydroxide solution, heat under a reflux condenser in a water bath for 30 minutes. Add 50 ml of 1 M sulphuric acid. Cool and filter, rinse the flask with 1 M sulphuric acid and combine the filtrate and washings in a 100-ml flask. Make up the volume with 1 M sulphuric acid. Label this solution as Solution A. Into two matched test tubes, introduce respectively, 1 ml of solution A and 1 ml of glycerin solution in 1 M sulphuric acid (Solution B). Add 0.5 ml of solution of sodium periodate to each tube, mix and allow to stand for 5 minutes. Add to each tube 1 ml of decolourized fuchsin solution and mix so that any precipitate disappears. Place the tubes in a beaker containing water at 40°C and observe for 10 to 15 minutes.

**A-10.2.1** The material shall be taken to have passed the test if any bluish violet colour in the tube containing Solution A is not more intense than that in the tube containing Solution B.

(PCD 19)